

Diameter(mm)	2,30	2,30	2,70	2,70	3,50	3,50	4,00	4,00	4,50	4,50
Thread Type	cortical	conical	cortical	cortical	cortical	cortical	cancellous	cancellous	cortical	cortical
Head Type	locking	locking	non-locking	locking	non-locking	locking	non-locking	locking	non-locking	locking
Hole Diameter	0	0	0	0	0	0	0	0	0	0
Thread Type	full	full	full	full	full	full	full	full	full	full
Drive Type	T7 Torx	T8 Torx	T15 Torx	T15 Torx	T15 Torx	T15 Torx	T15 Torx	T15 Torx	3,5 Hex	3,5 Hex
Total Screw Length Range (mm)	8-30	8-50	12-40	12-40	12-80	12-80	10-100	10-100	14-70	14-70
Diameter(mm)	6,50	6,50	6,50	6,50	6,50	6,50	4,00	4,00	4,00	6,50
Thread Type	cancellous	cancellous	cancellous	cancellous	cancellous	cancellous	cannulated cancellous	cannulated cancellous	cannulated cancellous	cannulated cancellous
Head Type	locking	locking	locking	non-locking	non-locking	non-locking	non-locking	non-locking	non-locking	locking
Hole Diameter	0	0	0	0	0	0	1,3 mm	1,3 mm	1,3 mm	2,6 mm
Thread Type	16 mm	32 mm	full	16 mm	32 mm	full	16 mm	32 mm	full	16 mm
Drive Type	3,5 Hex	3,5 Hex	3,5 Hex	3,5 Hex	3,5 Hex	3,5 Hex	2,5 Hex	2,5 Hex	2,5 Hex	3,5 Hex
Total Screw Length Range (mm)	30-90	45-90	30-90	30-90	45-90	30-90	10-70	30-70	20-70	30-90
Diameter(mm)	6,50	6,50	6,50	6,50	6,50	4,00	4,00	2,40	3,40	4,10
Thread Type	cannulated cancellous	cannulated cancellous	cannulated cancellous	cannulated cancellous	cannulated cancellous	malleolar	cannulated cortical	cannulated compression	cannulated compression	cannulated compression
Head Type	locking	locking	non-locking	non-locking	non-locking	non-locking	locking	locking	locking	locking
Hole Diameter	2,6 mm	2,6 mm	2,6 mm	2,6 mm	2,6 mm	0	1,3 mm	0,9	1,3	1,6
Thread Type	32 mm	full	16 mm	32 mm	full	32 mm	full	full	full	full
Drive Type	3,5 Hex	3,5 Hex	3,5 Hex	3,5 Hex	3,5 Hex	3,5 Hex	T15 Torx	1,5 Hex	2 Hex	2,5 Hex
Total Screw Length Range (mm)	30-90	30-90	30-120	50-120	30-120	30-70	30-95	8-30	16-30	16-50